

- 2014-01-09



B·R·A·H·M·S GmbH
Neuendorfstr. 25
16761 Hennigsdorf
GERMANY

Bulletin Number: **KB-169G** **Approval:** **jean-pierre bonneton - 2014-10-30** **Category:**
Instrument
Product Group: **KRYPTOR**
Model: All
Part Number:
Subject: KIM (Laboratory Information System)
Action: Information

Overview

THIS KRYPTOR BULLETIN REPLACES KB-169F WHICH IS NOT VALID ANYMORE

We added :

- a new warning related to Patient ID in the LIS provider documentation to explain that it is not possible to use the same patient ID for 2 different patient names.
- a new restriction for batch mode regarding the number of orders in the same frame.

We replaced the installation procedures for the different KRYPTOR systems (addition of a warning when switching from Pilote to KIM and modification of installation steps for KRYPTOR running under Windows 7).

The purpose of this Kryptor Bulletin is to guide you during a KIM installation and configuration. You can find different folders :

-[LIS provider documentation](#) : to forward to the LIS provider in order to do the configuration. This folder contains KIM Specifications and Host Lab instructions.



LIS provider documentation.zip

-[Installation Procedures for KIM software](#) : to guide you how to install KIM depending on the operating system running on KRYPTOR XPC.

2 sub-folders will be available :

- the first one with procedures for KRYPTOR compact / KRYPTOR compact PLUS running under Windows XP or Windows 7 and KRYPTOR running under Windows 98
- the second one for KRYPTOR running under Windows 7

For each subfolder, KIM installation procedure and Multi-sites specific procedures will be available.



KIM installation procedures.zip

-[License order](#) : to request to our spare part department the License that is mandatory to use the KIM software (KIM LICENSE ORDERING PROCEDURE.pdf and KIM License Number Request Form.pdf)



License order.zip

-[Procedures for KIM use and basic troubleshooting](#) : to explain the existing maintenance for KIM (LIS Maintenance.pdf) and give some advice how to resend result and troubleshoot if communication is interrupted (KIM Test Functions.pdf)



KIM use and troubleshooting.zip

-[KLIS 2.0.4](#) : to copy and paste on your XPC in order to install KIM software



KLIS 2.0.4.zip

Please find here below the confirmation form you have to fill in to confirm receipt of this Kryptor bulletin.

This confirmation form should be returned to the KRYPTOR Customer Service.



KB-169G_KIM (Laboratory Information System)_confirmation form.doc

Do not hesitate to contact Cezanne product support if you need further assistance.

KRYPTOR Customer Service,
Clinical Diagnostics Division
Thermo Fisher Scientific

280 Allée Graham Bell
30035 Nîmes cedex 1
Hotline France: +33 4 66 36 52 37
International Hotline: +33 4 66 36 52 46
Fax : +33 4 66 36 52 01

productsupport.brahms.frnim@thermofisher.com
www.thermoscientific.com

CEZANNE SAS, 280 Allée Graham Bell, Parc Scientifique Georges Besse, 30035 Nîmes, SAS
au capital de 1 120 000€, SIRET 432 032 712 00035