# Afinion™ Data Connectivity Converter

REF 1115660

User Manual	(US)	GB	2...9
Brukermanual	(NO)	10...17	
Bruksanvisning	(SE)	18...25	
Brugermanual	(DK)	26...33	
Käyttöohje	(FI)	34...41	
Bedienungsanleitung	(DE)	42...49	
Manuel d'utilisation	(FR)	50...57	
Guida per l'utente	(IT)	58...65	
Manual del usuario	(ES)	66...73	
Manual do Utilizador	(PT)	74...81	
Gebruikershandleiding	(NL)	82...89	
Εγχειρίδιο Χρήσης	(GR)	90...97	

**afinion**<sup>TM</sup>  
POINT OF CARE CENTER

**User Manual Conformity to Directive**

Afinion™ AS100 Analyzer equipped with Afinion™ Data Connectivity Converter (ADCC) meets all provisions in the European Directive 98/79/EC on *in vitro* diagnostic medical devices, and the ADCC component is CE marked accordingly.

**Safety standards**

Afinion™ AS100 Analyzer equipped with Afinion™ Data Connectivity Converter (ADCC) has been tested and found to be in conformity with IEC, UL, CAN/CSA-C22.2: 61010-1 (Safety requirements for electrical equipment for measurement, control, and laboratory use), IEC 61010-2-081:2001 + A1 and IEC 61010-2-101:2002 (Particular requirements for *in vitro* diagnostic (IVD) medical equipment).

**EMC standards**

Afinion™ AS100 Analyzer equipped with Afinion™ Data Connectivity Converter (ADCC) has been tested and found to be in conformity with EN 61326-1:2006 (Electrical equipment for measurement, control, and laboratory use – EMC requirements), EN 61326-2-6:2006 (*In vitro* diagnostic (IVD) medical equipment) and CFR 47: Telecommunications, Chapter I- FCC Part 15 – Radio Frequency Devices – Subpart B: unintentional radiators.

**TABLE OF CONTENTS**

1	Symbols and abbreviations	2
2	Intended use of the Afinion™ Data Connectivity Converter	3
3	Examining the package contents	3
4	Description of the Afinion™ Data Connectivity Converter	3
5	Installing the Afinion™ Data Connectivity Converter	4
5.1	Configuration of the ADCC	4
5.2	ADCC settings	4
5.3	Connecting ADCC	4
6	Web interface errors	6
7	Web interface ADCC software update	6
8	ADCC troubleshooting	7
8.1	LED status indicator signals	7
8.2	Confirm preparations	8
9	Testing procedures	8
10	Patient discretion	8
11	ADCC communication protocol data sheets	8
12	Maintenance and warranty	9
12.1	Cleaning and maintenance	9
12.2	Warranty	9
13	Disposal of the ADCC	9
14	Technical specifications	9
15	Third party software component	9

**1 SYMBOLS AND ABBREVIATIONS**

The following symbols and abbreviations are used in the product labelling and instructions for the Afinion™ Data Connectivity Converter. Please consult the Afinion™ AS100 Analyzer User Manual and Package Insert that comes with each Afinion™ test kit for symbols and abbreviations for the Afinion™ AS100 Analyzer system.

ADCC	Afinion™ Data Connectivity Converter
SN	Serial number
LOT	Lot number
LED	Light Emitting Diode
HIS	Hospital Information System
LIS	Laboratory Information System
ASTM	American Society for Testing and Materials
HL7	Health Level Seven
EPR	Electronic Patient Record
DHCP	Dynamic Host Configuration Protocol

IP	Internet Protocol
PAT	Patient
QC	Quality control
	CE marking of conformity to Directive 98/79/EC on <i>in vitro</i> diagnostic medical devices
	Manufacturer
	Storage temperature limitations
	Consult the ADCC User Manual
	Consult the Afinion™ AS100 Analyzer User Manual
	Waste Electrical and Electronic Equipment (WEEE)
	Power
	RS232
	Ethernet
	Warnings and precautions
	Universal Serial Bus (USB)

## 2 INTENDED USE OF THE AFINION™ DATA CONNECTIVITY CONVERTER

The ADCC is designed for Afinion™ AS100 Analyzer system customers at doctor's offices, hospitals or similar who need to connect the Analyzer to an electronic journal system (HIS/LIS).

The ADCC automatically transfers the Patient and Control results from the Afinion™ AS100 Analyzer to the HIS/LIS systems. The results are transferred to the HIS/LIS systems on the Ethernet in the standard formats HL7 or ASTM (selectable by configuration).

The ADCC is a component of the «Afinion™ AS100 Analyzer Connect», which comprises the Afinion™ AS100 Analyzer and the Afinion™ Data Connectivity Converter.

- The Afinion™ Data Connectivity Converter is only suitable for use with the Afinion™ AS100 Analyzer system.

## 3 EXAMINING THE PACKAGE CONTENTS

When unpacking, please check the contents against the list below and examine the components for signs of shipping damage.

The Afinion™ Data Connectivity Converter package includes:

- ADCC with RS232 and USB cable
- ADCC User Manual

If the package unit is found incomplete, please report missing items or shipping damage to your local supplier.

## 4 DESCRIPTION OF THE AFINION™ DATA CONNECTIVITY CONVERTER

The figure below shows the main exterior parts of the ADCC.

The ADCC may be reset to the default factory network settings by pushing a paperclip or pin into the Reset button.

## 5 INSTALLING THE AFINION™ DATA CONNECTIVITY CONVERTER

### 5.1 Configuration of the ADCC

Connect the USB cable to a USB port for power (i.e. the USB port on your PC).

Connect an Ethernet cable between your PC and the ADCC. To connect with the ADCC, the PC must have an Ethernet port in the same sub net as the ADCC. Manually select an IP address between 192.168.1.2 and 192.168.1.255. This can be selected in the control panel on your PC.

The ADCC web interface can be opened in a common web browser (i.e. Internet Explorer or Firefox) on your PC. Enter the ADCC IP address in the address bar from your web browser. The default IP address for the ADCC is 192.168.1.1.

The following screen will appear (with ADCC default settings):

Category	Parameter	Configuration
ADCC	DHCP	<input type="checkbox"/>
	IP address	192 [ ] 168 [ ] 1 [ ] 1 [ ]
	Network Mask	255 [ ] 255 [ ] 255 [ ] 0 [ ]
	Gateway	0 [ ] 0 [ ] 0 [ ] 0 [ ]
	HostName	ADCC12345
Remote System	IP address	0 [ ] 0 [ ] 0 [ ] 0 [ ]
	PortNo	0
Output	Communication	ASTM Low Level
	Results	New Results Only
	PAT / QC	Only Patient
Recipient	Receiving Application	EPR
	Receiving Facility	
	Send Patient ID as	HIS Patient ID
Password	Current	
	New	
	Repeat	
Reset	Configuration Settings to Default	<input type="checkbox"/>
	Clear Device List	<input type="checkbox"/>
		<input type="button" value="Update Settings"/>

Copyright ©Axis Shield 2009, All rights reserved.

Enter your preferred configuration settings (see 5.2 for an explanation of all settings).

Consult your network administrator and HIS/LIS administrator for required network settings.

When you have successfully updated all settings, click the «Update Settings» button in the web interface. The new settings will be stored in a non volatile memory and the ADCC box will automatically be restarted. The configured settings will then become operative.

### 5.2 ADCC Settings

Category	Parameter	Configuration
ADCC	DHCP enabled	<input checked="" type="checkbox"/> When using DHCP the IP address on the box will be assigned by the DHCP Server. No other network settings are necessary
	IP address	—
	Network Mask	—
	Gateway	—
	Host Name	Valid characters for hostname are [A-Z], [0-9], [-]. The Hostname length can be from 1-16 characters

Category	Parameter	Configuration
ADCC	DHCP disabled	<input type="checkbox"/> When using fixed IP address, the DHCP has to be disabled
	IP address	Please insert the fixed IP address [0-255/0-255/0-255/1-254]
	Network Mask	Please insert the network mask for the ADCC box [0-255/0-255/0-255/0-255]
	Gateway	Please insert the gateway [0-255/0-255/0-255/1-254]
	Host name	Valid characters for hostname are [A-Z], [0-9], [-]. The Hostname length can be from 1 -16 characters

Category	Parameter	Configuration
Remote System	IP address	Insert the IP address of the receiving system [0-255 0-255 0-255 1-254]
	Port No	[0-65535] (0 = not set)
Output	Communication	ASTM Low Level: For relevant information, see chapter 11. The communication protocol is based on ASTM E 1381 - 95 ASTM High Level: For relevant information, see chapter 11. The communication protocol is based on ASTM E 1394 - 97 HL7: For relevant information, see chapter 11. The communication protocol is based on HL7 version 2.4
	Results	New Results Only: Only test results obtained after first connection between ADCC and Afinion™ AS100 Analyzer will be transferred to the EPR system All Results: All results stored in the Afinion™ AS100 Analyzer memory and future test results obtained after first connection between ADCC and Afinion™ AS100 Analyzer will be transferred to the EPR system  If the ADCC has been connected to the Analyzer before, only results after last connection will be transferred. To send all results again the Clear device list must be used. See the category Reset
	PAT / QC	Only Patient: Only Patient results will be transferred to the EPR system Patient and QC: Patient and Control results will be transferred to the EPR system

Category	Parameter	Configuration
Recipient	Receiving Application	The length of name may consist of up to 30 characters. Invalid characters are \   & ^ ~ For relevant information, see chapter 11
	Receiving Facility (HL7 only)	The length of name may consist of up to 30 characters. Invalid characters are \   & ^ ~ For relevant information, see chapter 11
	Send Patient-ID as (HL7 only)	HIS Patient ID: For relevant information, see chapter 11 Visit Number: For relevant information, see chapter 11
Password	Current	The length of the password may consist of up to 15 characters
	New	New password
	Repeat	Confirm new password
Reset	Configuration settings to default	<input checked="" type="checkbox"/> Restore factory settings. Password and Device list will not be cleared
	Clear Device list	<input checked="" type="checkbox"/> It is possible to switch the ADCC to another Afinion™ AS100 Analyzer without changing the settings in the box. Choosing the «Clear Device list» will delete this memory
		«Update settings» will confirm and store new settings. The settings will be stored in a non volatile memory

## 5.3 Connecting ADCC

### Preparation of the Afinion™ AS100 Analyzer

The Afinion™ AS100 Analyzer should have SW 5.00 or higher installed. The SW version installed on the Afinion™ AS100 Analyzer will appear in the upper left corner of the screen during the Automatic start-up procedure. Contact your local Afinion™ Supplier for latest the Afinion™ AS100 Analyzer SW version. 

Connect the ADCC unit to the Afinion™ AS100 Analyzer serial port with the ADCC RS232 cable. Connect the ADCC Ethernet port to the HIS/LIS system network with an Ethernet cable. Power the ADCC by connecting the ADCCs USB cable to the Afinion™ AS100 Analyzer USB port and switch ON the Afinion™ AS100 Analyzer. The ADCC will now automatically transfer the results from the Afinion™ AS100 Analyzer to the HIS/LIS systems

 Allow the Afinion™ AS100 Analyzer to perform the Automatic start-up procedure. 

### Adjusting time and date

 Precaution for ADCC with SW 1.00 installed: The correct time and date should be set on the Afinion™ AS100 Analyzer before connecting to the ADCC. Later time adjustments should be done prior to running any tests that day. When returning to standard time from daylight saving time, allow a minimum of one hour from the last test run before adjusting the analyzer time. 

## 6 WEB INTERFACE ERRORS

Error text	Action to take	
Invalid IP configuration	IP address	[0-255/0-255/0-255/1-254]
	Network Mask	[0-255/0-255/0-255/0-255]
	Gateway	[0-255/0-255/0-255/1-254]
The password is incorrect	Password is missing or not correct	
The hostname contain invalid characters	Allowed characters are [A-Z], [0-9] and [-]	
The fields «Receiving Application» or «Receiving Facility» contains invalid characters	The field contains invalid characters \   & ^ ~	
The passwords entered are not the same. Please re-enter.	Please re-enter password	

## 7 WEB INTERFACE ADCC SOFTWARE UPDATE

Current SW version, User Manual and data sheets may be downloaded at Afinion.net.

Store the SW file locally on your computer. Open your ADCC web interface. The SW version is displayed in the upper right corner of the interface. Choose «ADCC Software Update».

Browse to the SW file. If your settings are password protected, type in your password.  
Click on «Start update».

The SW will transfer to the ADCC.

## Afinion™ Data Connectivity Converter

- ADCC Configuration
- ADCC Software Update

Software update will be performed.  
After successful software update you will be forwarded automatically to the configuration page!  
The update needs about 20 seconds. Please wait!

Copyright © Axis Shield 2009, All rights reserved.

The SW will then be updated on the ADCC. The update takes about 25 seconds.

## Afinion™ Data Connectivity Converter

- ADCC Configuration
- ADCC Software Update

Category	Parameter	Configuration
ADCC	DHCP	<input type="checkbox"/>
	IP address	<input type="text"/> 192 <input type="text"/> 168 <input type="text"/> 1 <input type="text"/> 1
	Network Mask	<input type="text"/> 255 <input type="text"/> 255 <input type="text"/> 255 <input type="text"/> 0
	Gateway	<input type="text"/> 0 <input type="text"/> 0 <input type="text"/> 0 <input type="text"/> 0
	HostName	<input type="text"/> ADCC12345
Remote System	IP address	<input type="text"/> 0 <input type="text"/> 0 <input type="text"/> 0 <input type="text"/> 0

SN ADCC12345

2.00

The SW update is complete when you are forwarded automatically to the configuration page. The new SW version will be displayed in the upper right corner.

- If the version in the upper right corner has not changed, the SW update failed. Repeat the ADCC SW update procedure. If this update fails, please contact your local Afinion™ supplier.
- Power must not be disconnected during SW update. Loss of power during SW update might lead to ADCC failure.
- During SW update of ADCC, the unit should be disconnected from the Afinion™ AS100 Analyzer. Tests performed on an Afinion™ AS100 Analyzer connected to an ADCC during SW update might lead to loss of result transfer to HIS/LIS system.

## 8 ADCC TROUBLESHOOTING

### 8.1 LED status indicator signals

Symbols				
Situations	Power LED (Green)	RS232 LED (Yellow)	Ethernet LED (Yellow)	Warning LED (Red)
Power off	OFF	OFF	OFF	OFF
Power on	Steady	Any state	Any state	Any state
Power on, no HIS/LIS communication. Wait 1 min. Check ADCC configuration and Ethernet connection	Steady	Any state	OFF	Any state
Power on, no Afinion™ AS100 Analyzer communication Analyzer in Automatic start up procedure. Wait 1 min after Analyzer has performed Automatic start-up procedure	Steady	OFF	Any state	Any state
Afinion™ AS100 Analyzer connection established	Steady	Steady	Any state	Any state
HIS/LIS connection established	Steady	Any state	Steady	Any state
Connection succeeded, HIS/LIS and Afinion™ AS100 Analyzer connection established	Steady	Steady	Steady	OFF
Afinion™ AS100 Analyzer communication failed. Wait 1 min. Analyzer in self-test. Wait 1 min. after Automatic start-up procedure is performed	Steady	Flashing	Any State	Flashing
HIS/LIS communication fatal error	Steady	Any State	Flashing	Flashing
Unrecoverable error not related to port communication	Steady	OFF	OFF	Flashing
SW upgrade is finalizing	Steady	LED is lit in sequence during upgrade. SW upgrade is successfully performed when LED flashes three times simultaneously		
Factory reset button activated	Single flash	Single flash	Single flash	Single flash

If the results are not transferred but the LED lights are not displaying any of the errors above, please reboot ADCC. If the error persists, please contact your local Afinion™ supplier.

For US customers a toll free line for technical support is available

## 8.2 Confirm preparations

These steps are not required but it is a path to confirm that the preparations were successful and/or for troubleshooting.

### IP address on PCs

Enter «CMD» in RUN in ms-dos window. Then enter the command «ipconfig» into the ms-dos window. This lists all the Ethernet ports and the addresses that are associated to these ports. If there are multiple Ethernet ports, the ones not used to configure the ADCC should be deactivated.

Consult your PC user manual for this procedure. Please find below an example of a result with IP address e.g. 192.168.1.100 and Subnet Mask 255.255.255.0

Ethernet adapter Local Area Connection:

IP Address. .... : 192.168.1.100

Subnet Mask. .... : 255.255.255.0

### Confirm connection between ADCC and PC

Connect the ADCC to the PC with an Ethernet cable and power from the Afinion USB. Confirm the connection between ADCC and PC with the ms-dos command «ping ADCC IP address» e.g. ping 192.168.1.1

#### The output text will be:

C:\Documents and Settings\vrfr>ping 192.168.1.1

Pinging 192.168.1.1 with 32 bytes of data:

Reply from 192.168.1.1: bytes=32 time=8ms TTL=64

Reply from 192.168.1.1: bytes=32 time<1ms TTL=64

Reply from 192.168.1.1: bytes=32 time<1ms TTL=64

Reply from 192.168.1.1: bytes=32 time<1ms TTL=64

Ping statistics for 192.168.1.1:

packets: Sent = 4, Received = 4, Lost = 0 (0% loss),

Approximate round trip times in milli-seconds:

Minimum = 0ms, Maximum = 8ms, Average = 2ms

The ADCC and the PC are now ready for ADCC configuration.

## 9 TESTING PROCEDURES

 Consult the Afinion™ AS100 Analyzer User Manual for general Test Cartridge handling instructions and how to analyse controls and patient samples on the Afinion™ AS100 Analyzer. Consult the Package Insert that comes with each Afinion™ test kit for assay specific information.

 •When running a patient test on an Afinion™ AS100 Analyzer connected to a HIS/LIS system the Patient ID function should be enabled and used. 

## 10 PATIENT DISCRETION

Configuration of the ADCC allows password protection of settings. For relevant information, see chapter 5.

 •When you export data that contains patient information, it is your responsibility to comply with your local regulations on protection of personal health information.

## 11 CONNECTING THE ADCC TO A LIS/HIS SYSTEM

The Afinion™ Data Connectivity Converter Protocol Data Sheet ASTM and HL7 should be downloaded at [www.afinion.net](http://www.afinion.net) or contact your local Afinion™ supplier.

 •The selected Afinion™ Data Connectivity Converter Protocol Data Sheet must be examined carefully before implementing the ADCC.

 •Precaution for ADCC with SW 1.00 installed: Precaution for ACR. The ACR observation value shall be interpreted as not valid if albumin and/or creatinine are below or above measurement limits, independent of the content in ACR abnormal flag field. See section 4.6 in ASTM protocol SW 1.00 and 4.8 in HL7 protocol SW 1.00 for further details.

## 12 MAINTENANCE AND WARRANTY

### 12.1 Cleaning and maintenance

No maintenance of the Afinion™ Data Connectivity Converter is required other than cleaning the exterior. Clean the outside of the ADCC with clean, lint-free and non-abrasive cloth dampened in water or a mild detergent.

### 12.2 Warranty

Axis-Shield PoC guarantees a warranty period of 12 months from the date of delivery to the user. For this period it is guaranteed that the Afinion™ Data Connectivity Converter works within quality requirements valid on the date of manufacture release.

Axis-Shield PoC is prepared to compensate for faulty material or malfunctions within the limitations of the warranty. The warranty is limited to providing a replacement Afinion™ Data Connectivity Converter.

Tampering with the internal components, damage due to operating errors, misuse and overlooking essential information with respect to warnings and precautions described in this User Manual will invalidate the warranty.

The acknowledgment of claims shall immediately be reported to your local Afinion™ supplier.

## 13 DISPOSAL OF THE ADCC

For correct disposal according to the Directive 2002/96/EC (WEEE), contact your local Afinion™ supplier.

## 14 TECHNICAL SPECIFICATIONS

Afinion™ Data Connectivity Converter	
Size	110 x 248 x 372 mm
Weight	125 g
Power supply	The ADCC gets power (5 V) from the Afinion through the USB cable. The USB cable acts, in this case, only as a power connector
Communication interface	Afinion™ AS100 Analyzer to ADCC: RS232 ADCC output: Ethernet
Operation conditions	
Temperature	15 – 40 °C / 59 – 104 °F
Relative humidity	0 – 90%
Storage and transport (in the original container)	
Temperature	-40 – 70 °C / -40 – 158 °F
Relative humidity	0 – 90%

## 15 THIRD PARTY SOFTWARE COMPONENT

This product contains third party lwIP TCP IP Stack software.

<http://www.sics.se/~adam/lwip/>

Copyright (c) 2001, 2002 Swedish Institute of Computer Science.  
All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:  
1. Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.  
2. Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.  
3. The name of the author may not be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE AUTHOR «AS IS» AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE AUTHOR BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

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Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

1. Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
2. Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.
3. The name of the author may not be used to endorse or promote products derived from this software without specific prior written permission.

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Category	Parameter	Configuration
ADCC	DHCP	<input type="checkbox"/>
	IP address	192 [ ] 168 [ ] 1 [ ] 1 [ ]
	Network Mask	255 [ ] 255 [ ] 255 [ ] 0 [ ]
	Gateway	0 [ ] 0 [ ] 0 [ ] 0 [ ]
	HostName	ADCC12345
Remote System	IP address	0 [ ] 0 [ ] 0 [ ] 0 [ ]
	PortNo	0
Output	Communication	ASTM Low Level
	Results	New Results On
	PAT / QC	Only Patient
Recipient	Receiving Application	EPR
	Receiving Facility	
	Send PatientID as	HIS Patient ID
Password	Current	
	New	
	Repeat	
Reset	Configuration Settings to Default	<input type="checkbox"/>
	Clear Device List	<input type="checkbox"/>
		<b>Update Settings</b>

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ADCC	Afinion™ Data Connectivity Converter
SN	Seriennr.
LOT	Lotnr.
LED	Lysdiode
HIS	Hospital Information System (patientjournalsystem)
LIS	Laboratory Information System (laboratoriejournalsystem)
ASTM	American Society for Testing and Materials
HL7	Health Level Seven
EPR	Electronic Patient Record (elektronisk patientjournal)
DHCP	Dynamic Host Configuration Protocol

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Category	Parameter	Configuration
ADCC	DHCP	<input checked="" type="checkbox"/>
	IP address	192 168 1 1
	Network Mask	255 255 255 0
	Gateway	0 0 0 0
	HostName	ADCC12345
Remote System	IP address	0 0 0 0
	PortNo	0
Output	Communication	ASTM Low Level
	Results	New Results Only
	PAT / QC	Only Patient
Recipient	Receiving Application	EPR
	Receiving Facility	
	Send PatientID as	HIS Patient ID
Password	Current	
	New	
	Repeat	
Reset	Configuration Settings to Default	<input type="checkbox"/>
	Clear Device List	<input type="checkbox"/>

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The screenshot shows the main configuration interface for an ADCC device. On the left, there's a sidebar with navigation links: 'ADCC Configuration' and 'ADCC Software Update'. The main area displays a table with columns 'Category', 'Parameter', and 'Configuration'. For the 'ADCC' category, the parameters listed are 'DHCP' (unchecked), 'IP address' (with values 192, 168, 1, 1), 'Network Mask' (with values 255, 255, 255, 0), 'Gateway' (with values 0, 0, 0, 0), and 'HostName' (set to ADCC12345). For the 'Remote System' category, the 'IP address' parameter is shown with values 0, 0, 0, 0. In the top right corner of the main window, there's a red circle highlighting the 'SN ADCC12345' and '1.00' text, which indicates the serial number and software version respectively.

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This screenshot shows the 'Software Update' screen. It features a sidebar with 'ADCC Configuration' and 'ADCC Software Update' links. The main part of the screen has fields for 'File:' (with a 'Browse...' button) and 'Password:' (with a password input field). Below these fields is a 'Start Update' button. At the bottom, a copyright notice reads 'Copyright ©Axis Shield 2009, All rights reserved.'

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This screenshot is similar to the previous one, showing the 'Software Update' screen. The 'File:' field now contains the path 'C:\Documents and Settings\gspe\Desktop\ADCC SW'. The 'Start Update' button is visible at the bottom, along with a progress indicator showing '10 %'. The copyright notice at the bottom is identical to the previous screenshot.

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Category	Parameter	Configuration
ADCC	DHCP	<input checked="" type="checkbox"/>
	IP address	192 [ ] 168 [ ] 1 [ ] 1 [ ]
	Network Mask	255 [ ] 255 [ ] 255 [ ] 0 [ ]
	Gateway	0 [ ] 0 [ ] 0 [ ] 0 [ ]
	HostName	ADCC12345
Remote System	IP address	0 [ ] 0 [ ] 0 [ ] 0 [ ]
	PortNo	0
Output	Communication	ASTM Low Level
	Results	New Results Only
	PAT / QC	Only Patient
Recipient	Receiving Application	EPR
	Receiving Facility	
	Send Patient ID as	HIS Patient ID
Password	Current	
	New	
	Repeat	
Reset	Configuration Settings to Default	<input type="checkbox"/>
	Clear Device List	<input type="checkbox"/>
		<b>Update Settings</b>

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Category	Parameter	Configuration
ADCC	DHCP	<input checked="" type="checkbox"/>
	IP address	192 [ ] 168 [ ] 1 [ ] 1 [ ]
	Network mask	255 [ ] 255 [ ] 255 [ ] 0 [ ]
	Gateway	0 [ ] 0 [ ] 0 [ ] 0 [ ]
	HostName	ADCC12345
Remote System	IP address	0 [ ] 0 [ ] 0 [ ] 0 [ ]
	Port No	0
Output	Communication	ASTM Low Level
	Results	New Results Only
	PAT / QC	Only Patient
Recipient	Receiving Application	EP R
	Receiving Facility	
	Send PatientID as	HIS Patient ID
Password	Current	
	New	
	Repeat	
Reset	Configuration Settings to Default	<input type="checkbox"/>
	Clear Device List	<input type="checkbox"/>
<input type="button" value="Update Settings"/>		

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The screenshot displays the configuration interface for the Afinion™ Data Connectivity Converter. The main window title is "Afinion™ Data Connectivity Converter". On the left, there's a sidebar with icons for "ADCC Configuration", "ADCC Software Update", and "Axis-Shield". The main panel contains several sections with input fields and dropdown menus:

- Category**: ADCC
- Parameter**: IP address, Network Mask, Gateway, Host Name, Remote System IP address, Port No., Output Communication (set to ASTM Low Level), Results (set to New Results Only), PAT / QC, Recipient Receiving Application (set to EP R), Recieving Facility, Send PatientID as (set to HIS Patient ID), Password Current, New, Repeat, Reset Configuration Settings to Default, Clear Device List.
- Configuration**: Includes fields for IP address (192.168.1.1), Network Mask (255.255.255.0), Gateway (0.0.0.0), Host Name (ADCC12345), and a MAC address field (ADCC12345).

At the bottom right of the interface, it says "Copyright ©Axis Shield 2009, All rights reserved."

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