


# Service Manual


***Rapidchem™*** 700

**Bayer** 



**REF** 04663851  
©2003 Bayer Corporation  
(all rights reserved)

Flow Chart



## Table of Symbols

---



Important information



Conditions which may cause data loss or analyzer malfunction



Biohazard warning



Electrical hazard

# Service Manual Contents

---

Software Flow Chart . . . . .	inside front cover
Table of Symbols . . . . .	i

## Product Description

System Overview . . . . .	1
Sensor Module . . . . .	2
Sampler . . . . .	2
Valve Module . . . . .	2
Printer Assembly . . . . .	3
CPU Assembly . . . . .	3
Pre-Amp PC Board . . . . .	3
Pump Assembly . . . . .	3
Contact Housing . . . . .	4
Display Assemblies . . . . .	4
Valve Module Cable . . . . .	5
Rear Housing Assembly . . . . .	5

## Troubleshooting

Troubleshooting Overview . . . . .	7
Troubleshooting Guide . . . . .	7
AC Power Verification . . . . .	8
DC Power Verification . . . . .	9
Display Assembly Verification . . . . .	10
Flash Device Verification . . . . .	11
Printer Verification . . . . .	12
Pump Verification . . . . .	13
Valve Verification . . . . .	13
Sensor Signal Input Verification . . . . .	14
Serial Communication Verification . . . . .	15

## Repair

Repair Overview . . . . .	16
Required Tools . . . . .	16
Rear Housing Removal . . . . .	17
Rear Housing Installation . . . . .	17
Printer Assembly Removal . . . . .	18
Printer Assembly Installation . . . . .	18
Display Assembly Removal . . . . .	19
Display Assembly Installation . . . . .	19
CPU Assembly Removal . . . . .	20
CPU Assembly Installation . . . . .	20
Flash Device Removal . . . . .	21
Flash Device Installation . . . . .	21
Printer EPROM Removal and Installation . . . . .	22
Access Door Removal . . . . .	23
Access Door Installation . . . . .	23
Pre-Amp PC Board/Contact Housing Removal . . . . .	24
Pre-Amp PC Board/Contact Housing Installation . . . . .	25
Valve Cable Removal . . . . .	26
Valve Cable Installation . . . . .	26
Pump Assembly Removal . . . . .	27
Pump Assembly Installation . . . . .	27
Sampler Removal . . . . .	28
Sampler Installation . . . . .	28
Software Upgrade Kit . . . . .	29
Appendixes: CPU Board Schematics . . . . .	A-J
Contacting Bayer Diagnostics . . . . .	K-N

## System Overview

The repairs and procedures in this manual are intended for authorized field service engineers. Instructions describe the repair or replacement of all modular internal components. It is recommended that the field service engineer also refer to the Operator's Manual for detailed instructions on operating and troubleshooting the analyzer to determine if repair or replacement is necessary for an internal component. The Replacement List indicates assemblies and components available to the field service engineer. Refer to the Repair section of this manual for the CPU schematics.



*All components should be considered biohazardous (contaminated with HIV or other pathogens).*

*Any replaceable component which comes in contact with biological samples, including the sampler, sensors, sensor module, pump tubing, valve module and reagent module may contain contaminated material. Treat all components, during use and disposal, as you would any biohazardous material.*





## Sensor Module

Contains sampler position sensors, bubble detector, door sensor, sampler, and compression plate. The Sensors  $\text{Na}^+$ ,  $\text{K}^+$ ,  $\text{Cl}^-$  or  $\text{Li}^+$ , spacer, and reference sensors are housed in the sensor module.




## Sampler


Attached to sensor module. Contains sample probe with probe wiper. The sampler is removable from the sensor module as a replacement part.



## Valve Module


Controls reagent (from reagent module) or air selection. The valve module contains a ceramic valve, motor, position detector. There are no serviceable parts inside the valve module.






### Printer Assembly

Contains the printer mechanism, printer software, and printer electronics. Prints analysis, calibration and diagnostic data.




### CPU Assembly

Contains transformer, software, inputs (sensors), motor control circuitry, and 2 serial ports.



### Pre-amp PC Board

Connects contact housing to the CPU assembly and amplifies sensor signals.



### Pump Assembly

Contains peristaltic pump with wire connector.



## Contact Housing


Contains all contacts from the sensors, as well as the contact housing cable for connection to the CPU assembly.




## Display Assemblies


(Analyzer operator interface.) Contains operator keypad and LCD display.

Numeric display




Graphic display  
(required for certain languages)





### Valve Module Cable

Connects valve module to the CPU board.



### Rear Housing Assembly

Rear housing contains the printer paper bay and cover, serial number label, and the openings for the fuse holder, AC power cord socket, battery connection, RS-232 serial interface connection, barcode reader connection, and power voltage setting.

## Troubleshooting Overview

This section describes troubleshooting the mechanical and electrical functions of the Rapidchem analyzer. Bayer recommends that dealers maintain an inventory of operational test assemblies for troubleshooting purposes. These test assemblies can substitute the questionable assemblies to confirm failures. For information on troubleshooting the fluidic system and sensors, refer to the Rapidchem Operator's manual.

- ▲ Run quality control material after any troubleshooting procedure.
- ▲ When removing or handling PC boards, it is necessary to be properly grounded to prevent static damage.

## Troubleshooting Guide

Determine the problem and corrective action needed.


Problem/Description	Problem Source	Corrective Action Steps
Blank display/No response	Power malfunction Defective display	AC power verification DC power verification Display assembly verification Flash Device verification
Printer does not print	Printer malfunction	Printer verification
No pump movement	Pump motor malfunction	Pump verification
No valve movement	Valve motor malfunction	Valve verification
Sensor voltages out of specification	Sensors or signal input malfunction	Check sensor module installation and performance Sensor signal input verification
Serial transmission failure	Serial cable Transmission Protocol	Serial communication verification

## AC Power Verification

When the power is first turned on, the Rapidchem analyzer will “beep,” and the display will light up. If this fails to happen, follow the verification procedure in this section. Review this entire section before replacing any components. First confirm that the power cord is inserted securely inside the AC power cord socket.






**Warning: Careless handling of AC power can be hazardous or fatal. Only qualified personnel should work with AC power.**



Factory-set Power Voltage Setting

- 1 Remove the power cord from the AC power cord socket.
- 2 Verify that the voltage selector switch voltage setting (110V or 220V) matches the power outlet voltage. If it does not, correct the setting.
- ★ 3 Verify that the fuses have the correct rating. Push down on the fuse holder release tab and pull the fuse holder straight out. Inspect both fuses.
  - For 110V analyzers, the fuses must be rated 0.8A, 250 volts
  - For 220V analyzers, the fuses must be rated 0.4A, 250 volts





Remove the fuses from the fuse holder. Using a standard ohmmeter, measure continuity between the ends of each fuse. Replace defective fuses.

If the voltage setting and the fuses are correct, but "no power" or "blank display" continues, proceed as follows:

- 4 Remove the rear housing. See *Rear Housing Removal* in the **Repair** section for instructions.
- 5 Reinstall the power cord (cord not shown).
- 6 Measure the voltage between the bottom sides of components L7 and L8 on the CPU assembly using a standard voltmeter. Refer to picture.
  - For 110V analyzers, the measured voltage must be  $100/120V \pm 10\%$
  - For 220V analyzers, the measured voltage must be  $230V \pm 10\%$

If no voltage, check the power cord and outlet. If the measured voltage is incorrect, replace the CPU assembly. See *CPU Assembly Removal / Installation* in the **Repair** section for instructions.

#### Test Point Voltage

Table 1

X1	+ 5.00V $\pm$ 0.15V
X4	+ 24.00V $\pm$ 5.00V
X5	+ 19.00V $\pm$ 4.00V
X6	- 19.00V $\pm$ 4.00V
X7	+ 5.00V $\pm$ 0.25V
X8	GND Reference
X9	- 12.00V $\pm$ 0.25V
X10	+ 12.00V $\pm$ 0.25V

Table 2

X2	GND Reference
X3	+ 8.5V $\pm$ 2.00V

Table 3

TP9	GND Reference
TP10	+ 2.000V $\pm$ .003V
TP12	- 12.00V $\pm$ 0.25V
TP13	- 5.00V $\pm$ 0.25V
TP14	+ 5.00 $\pm$ 0.25V
TP15	+ 15.00V $\pm$ 1.00V

#### DC Power Verification

- 1 Remove the rear housing. See *Rear Housing Removal* in the **Repair** section for instructions.
  - 2 Remove sensor module.
  - 3 Connect the power cord to an electrical outlet, then to the AC power cord socket.
  - 4 Connect the reference side of the voltmeter to X8 on the CPU board.
  - 5 Verify the voltage at each test point listed in Table 1.
  - 6 Connect the reference side of the voltmeter to X2 on the CPU board.
  - 7 Verify the voltage at each test point listed in Table 2.
  - 8 Connect the reference side of the voltmeter to TP9 on the pre-amp PC board.
  - 9 Verify the voltage at each test point listed in Table 3.
- ★ No adjustments are permitted. Failed boards should be replaced. See *CPU Assembly Removal / Installation* in the **Repair** section for instructions.

## Display Assembly Verification


For blank display, or faulty keypad operation, proceed as follows:

- 1** Remove the power cord from the AC power cord socket.
- 2** Remove the rear housing. *See Rear Housing Removal in the **Repair** section for instructions.*
- 3** Disconnect and reinstall the display cable to the CPU P4 connector. *See Display Assembly Removal / Installation in the **Repair** section for instructions.*
- 4** Reinstall the power cord. If the display now functions, remove the power cord, install the rear housing, reinstall the power cord, and resume normal operation. *See Rear Housing Installation in the **Repair** section for instructions.* If the display or keypad remain faulty, proceed with step 5.
- 5** Remove the power cord from the AC power cord socket. Connect a good "test" display assembly to the CPU assembly. Reinstall the power cord.

If the "test" display functions properly, the original display is defective. Replace with a new display assembly.

If the "test" display does not function, the CPU is defective. Install a new CPU assembly.


*See CPU Assembly Removal / Installation in the **Repair** section for instructions.*



## Flash Device Verification

If the Flash Device is not completely inserted into position U3 on the CPU, the Rapidchem display assembly will remain blank on power up, and the analyzer will appear "off". Inspect the insertion. Verify that it is properly installed by pressing on the device with your thumb.

If it is necessary to inspect the Flash Device itself, remove the Flash Device with the removal tool by prying each corner of the Flash Device until it can be removed by hand. Inspect the Flash Device's metal connectors. They should appear curled and separated from each other. If any of the connectors appear bent, gently straighten them. If the Flash Device is permanently damaged, contact your Rapidchem dealer.



Reinstall the Flash Device into the U3 position with the chamfered edge in the lower right corner. Incorrect Flash Device orientation may damage the Flash Device and/or socket.




## Printer Verification

If the printer is not actively printing:


▲ Verify that the printer paper is correctly installed.

- 1 Proceed to the **SETUP MENU? USER OPTIONS?** is displayed. Press NO. Press YES to **CONFIGURATION?** and **PRINTER STATUS?**. Confirm that the printer is ON by pressing YES or NO to the displayed message. Confirm that the paper is threaded properly.
- 2 Press 9 on the keypad. The printer head should move across the paper. If the printer head moves, resume normal operation. If the printer head does not move, proceed to Step 3.
- 3 Remove the power cord from the AC power cord socket, wait 10 seconds, then reinstall the power cord. If the printer is still not functioning, remove the power cord and the rear housing, then remove the printer assembly. Inspect for paper jams or obstructions in the printer head area. See *Rear Housing Removal and Printer Assembly Removal* in the **Repair** section for instructions.
- 4 Reinstall the printer assembly. See *Printer Assembly Installation* in the **Repair** section for instructions. Confirm that the printer cable is completely inserted into the P3 connector on the CPU. Reinstall the power cord. If the printer fails to print, install a known working "test" printer assembly in the analyzer.



If the test printer assembly works, install a new printer assembly in the analyzer.

If the test printer assembly fails, the CPU is defective. Install a new CPU assembly. See *CPU Assembly Removal / Installation* in the **Repair** section for instructions.




## Pump Verification



For suspected pump problems, proceed to **PUMP TEST?** under **TEST COMPONENTS?** in the **DIAGNOSTICS?** section of the **SECOND MENU**. Press YES. The screen displays **ROTATING . . .**. Look for pump rotation. If full pump rotation occurs, **FULL ROTATION?** is displayed. Press YES and **PASS** is displayed.

When there is no pump movement or partial pump movement, press NO and **FAIL** is displayed. Remove the pump tubing and repeat the **PUMP TEST?**. If the test passes, install a new pump tubing and resume operation. If the test fails, install a new pump assembly. See *Pump Assembly Removal / Installation* in the **Repair** section for instructions.



## Valve Verification


For suspected valve problems, proceed to the **VALVE TEST?** under **TEST COMPONENTS?** in the **DIAGNOSTICS?** section of the **SECOND MENU?** Press YES. The screen displays **ROTATING . . .**, then **FULL ROTATION?**. Look for rotation of the valve knob and press YES or NO. When full valve knob rotation occurs, press YES and **TESTING . . .** is displayed. Internal testing verifies correct valve positioning for the Cal A, Cal B, Rinse, Air and Home positions. **PASS** or **FAIL** is displayed.

When there is no valve knob movement or partial valve knob movement, press NO and **FAIL** is displayed. Repeat the **VALVE TEST?** If the test fails again, install a new valve module. See **Rapidchem Operator's Manual** for valve module replacement instructions.

## Sensor Signal Input Verification

Rapidchem sensors that report out of range voltages may be defective and require replacement. Sensor troubleshooting is discussed in the Rapidchem Operator's Manual. If the Rapidchem analyzer reports repeated sensor voltage errors, follow this procedure.

- 1 Remove the Rapidchem sensors from the sensor module and check for any fluid leakage around the sensors and their pins. Examine the contact housing connectors for any corrosion. All of these surfaces must be completely dry and salt free. Reinstall the sensors and attempt recalibration. If the voltage problems continue, it is likely that a defect exists in the amplifier section of the CPU assembly, or the contact housing is defective.
- 2 Remove the rear housing. See *Rear Housing Removal* in the **Repair** section for instructions.
- 3 Install a new pre-amp PC board. See *Pre-amp PC Board/Contact Housing Removal / Installation* in the **Repair** section for instructions.



## Serial Communication Verification

- 1 Using a Rapidchem serial cable, connect the analyzer to a device with a known working serial port.
- 2 Refer to the serial interface specifications on this page.
- 3 When correctly linked, the data that is on the external device should match the data on the printout.

▲ See Rapidchem Operator's Manual for Computer Connection setup.

▲ The data is always sent to the port on the rear of the analyzer at a baud rate of 2400 bps. The bit configuration is 8 data bits, 1 stop bit, no parity. If the data is not being transmitted, the CPU board will need to be replaced.

Serial Interface Specifications	
Connector	25 or 9 pin or a D-sub male or both
Pin Type	RS-232
Pin Assignments:	
Pin 1	Chassis, Not connected
Pin 2	TxD
Pin 3	RxD
Pin 4	RTS
Pin 5	CTS
Pin 7	Signal Ground
Transmission Rate	2400 bits per second
Bit Configuration	8 data bits, 1 stop
Bit Parity	None
Control Codes	None
All lines end with carriage return line feed (CRLF)	
All transmission is ASCII string	


## Repair Overview

To replace Rapidchem assemblies, disassemble the analyzer in the order described in this section, referring to the photographs for assistance. After rear housing removal, the printer assembly and display assembly can be removed without removing the CPU. Access to all other assemblies requires removing the CPU board.

Before performing major repairs, proceed to **PRINT SETUP?** under **STORED INFORMATION?** in the **SECOND MENU?** and print the following settings:

**USER OPTIONS?,  
CONFIGURATION?,  
PATIENT LIMITS?,  
PATIENT INFORMATION?**


▲ Quality control data and patient results will be lost if the CPU board is changed.



### Required Tools

Only 4 tools are required to remove or install the Rapidchem electronic and mechanical assemblies:


- Phillips screwdriver, #2 magnetic
- Phillips screwdriver, #1
- Needle nose pliers
- Flash Device Removal Tool (provided with software replacements and updates)



## Rear Housing Removal

Cut the printer paper to separate it from the analyzer and remove paper. Remove the 5 screws from the rear housing.

Pull back the rear housing.



## Rear Housing Installation


Align the rear housing with the front housing. Verify that all cables are enclosed inside the housings. Insert and tighten the 5 screws into the rear housing. Do not over tighten.

Reinstall the printer paper.




## Printer Assembly Removal

Disconnect the P1 connector from the CPU assembly.




Slide the printer assembly out of the housing tracks. Avoid pinching the flex cable under the housing.



## Printer Assembly Installation

Slide the printer assembly into the front housing cavity, using a thin, flat object (e.g., business card) to guide the printer flex cable under the housing. Avoid pinching the cable.

Insert the printer cable into the P1 connector on the CPU as shown above.







## Display Assembly Removal

The reagent module, sensor module, valve module, and printer assembly must be removed before removing the display assembly.


Disconnect the P2 connector from the CPU assembly.



Remove the 4 screws that run along the top of the rear surface of the front housing.



After screw removal, grasp the display assembly as shown and gently lift up and away from the front housing, guiding the ribbon cable through the opening in the front housing.



## Display Assembly Installation

Insert the display cable through the rectangular window in the front housing.

Align the display assembly insertion hooks with the openings in the front housing. Push the assembly against the front housing until it locks. Enclose the top access door post under the upper right corner of the display housing.

Reinstall 4 screws along the top surface of the front housing. Attach the display cable to the P2 connector on the CPU assembly.





## CPU Assembly Removal

★ Before replacing the CPU assembly, you must have a new CPU assembly and a Flash Prom Replacement Kit. It is necessary to go through the steps listed in the next section, "Flash Device Removal." Prior to removing the old CPU board, install the new CPU board with the new Flash Prom and follow the steps listed in the next section, "Flash Prom Installation."

Remove the printer assembly before removing CPU assembly. Disconnect cables at CPU positions P1, P2, P5, P7, P6 and P13.


Remove the screw in the center of the CPU (P15), while supporting the transformer. Carefully slide the CPU assembly off the 4 corner mounting posts and the pre-amp PC board.



## CPU Assembly Installation

Carefully line up the four corner mounting holes on the CPU assembly with the mounting posts on the front housing. Gently push in the CPU assembly, making the connection between the pre-amp PC board and CPU assembly.

Secure CPU assembly to front housing with mounting screw at P15.




## Flash Device Removal



Before changing the flash device on the Rapidchem, follow these steps:

(Failure to follow these steps will result in analyzer malfunction.)

STANDBY?	Press yes
STANDBY IN 1 MIN.	Press yes, yes, yes, 9, yes, yes
SERVICE MENU?	Press yes
PROG. DOWNLOAD?	Press no
LOADER DOWNLOAD?	Press no
REPLACE PRE-AMP?	Press no
REPLACE FLASH?	Press yes
PLEASE CONFIRM!	Press yes
POWER DOWN NOW	Disconnect power



Remove the flash device with the removal tool by prying each corner of the flash device until it can be removed by hand.

## Flash Device Installation

Install the flash device into the U3 position with the chamfered edge in the lower right corner. Incorrect flash device orientation may damage the flash device and/or socket.



After changing the flash device on the Rapidchem, power up and follow these steps:

STANDBY? Press yes


STANDBY IN 1 MIN. Press yes, yes, yes, yes, yes

QC MENU? Press yes

UNIT SETTINGS? Press no

CLEAR MEMORY? Press yes

PLEASE CONFIRM! Press yes



### Access Door Removal

Remove the CPU assembly and loosen the top left screw on the rear of the analyzer.



Lift the right corner of the display assembly to disengage the access door post. Lift the door to remove the bottom post.

### Access Door Installation

Insert the bottom door post into the front housing base. Lift the right corner of the display assembly and insert the access door post. Tighten the screw on the rear of the front housing.

- ▲ In order to properly transfer data stored from the existing pre-amp PCB to a new pre-amp PCB, it is necessary to follow the steps listed in the chart below. (Failure to follow these steps will result in a “MISMATCH ERROR” on the printout.)

★ Before changing the Pre-amp PC board on the Rapidchem, follow these steps:



STANDBY?	Press yes
STANDBY IN 1 MIN.	Press yes, yes, yes, 9, yes, yes
SERVICE MENU?	Press yes
PROG. DOWNLOAD?	Press no
LOADER DOWNLOAD?	Press no
REPLACE PRE-AMP?	Press yes
PLEASE CONFIRM!	Press yes
POWER DOWN NOW	Disconnect power


### Pre-amp PC Board/Contact Housing Removal

Remove the screw which connects the grounding strap to the pre-amp PC board.

Remove the 6 screws from the pre-amp PC board.

Remove the pre-amp PC board with the contact housing attached from the front housing. The contact housing can be removed from the pre-amp PC board by gripping both sides and pulling apart.





### Pre-amp PC Board/Contact Housing Installation

Align the opening on the pre-amp PC board with the CPU board mounting post. This will align the screw holes. Insert the 6 screws to secure. Do not overtighten. Connect the grounding strap to the pre-amp PC board.



After changing the Pre-amp PC board on the Rapidchem power up and follow these steps:

STANDBY?	Press yes
STANDBY IN 1 MIN.	Press yes, yes, yes, yes, yes
QC MENU?	Press yes
UNIT SETTINGS?	Press no
CLEAR MEMORY?	Press yes
PLEASE CONFIRM!	Press yes


## Valve Cable Removal

Remove the retaining clips securing the valve cable with needle nose pliers.

Caution: DO NOT damage the posts.

By squeezing the retaining clips with the pliers, the clips should lift off the posts.

Remove the valve cable from the front housing.




## Valve Cable Installation

Insert the valve cable into the back side of the front housing so that the ribbon cable is positioned on the left side of the connector.

Install 2 new retaining clips (top and bottom) on the cable posts to secure the connector. Push clips on using needle nose pliers.






## Pump Assembly Removal



Remove the CPU assembly before removing the pump assembly. Using a #1 Phillips, remove the 2 screws and washers from the front side of the front housing.



Remove the screws which connect the grounding strap to the pre-amp PC board and the pump assembly. Place the grounding strap aside. Slide the pump assembly through the back side of the front housing.



## Pump Assembly Installation

Insert the pump through the backside of the front housing. Align the pump bracket holes with the 2 mounting posts on the front housing so that the pump motor cable is positioned below the pump assembly.

While holding the pump assembly in position, insert the 2 screws with washers through the front side of the front housing to secure the pump in position.

Attach the grounding strap to the pump assembly and the pre-amp PC board.




## Sampler Removal




If the sampler requires replacement, it can be removed from the sensor module. Using a coin or fingers, remove the retaining screw from the sensor module.

- ▲ Also refer to package insert for sampler removal and installation updates.



Remove the screw mount (includes spring).

Rotate the sampler into the capillary position, then pull it out of the sensor module.




## Sampler Installation

Install the new sampler (includes sample probe and wiper) in the sensor module. Insert the sampler alignment pin into the sensor module. Rotate the sampler into the closed position.

Install the screw mount (includes spring) into the sampler. Replace the retaining screw and tighten.

- ▲ Do not overtighten the retaining screw.



## Software Upgrade Kit

This kit includes a 3.5" Floppy Disk and a serial cable.

## Analyzer Setup

Using the serial cable, connect to the PC serial port. Connect the RJ11 plug/connector to the barcode port on the rear panel of the Rapidchem.

Store files provided by Bayer (\*.hex and \*.key) in a convenient directory. The name of the hex file is the CRC number used to validate the successful update of the program.

## PC Setup

Open Hyper Terminal under Programs > Accessories > Communications.

Enter a new connection name and click OK. Select Direct to Com port connected and click OK.

Change to 4800 bits per second, 8 Data Bits, Parity NONE, Stop Bit 1, and Flow Control Xon/Xoff.

Click at Advanced, uncheck Use FIFO Buffer, and press OK twice.

Click at File, then select Properties to confirm all settings that selections above are correct.

Select File, then Save As, and enter a new file name.

Load New User Program

Open Hyper Terminal with your new hyper terminal file.

Turn analyzer on. Select **STANDBY?** under **SECOND MENU?**

When the **STANDBY IN 1 MIN.** message is blinking, press yes three times, press 9 once, then press yes twice.

**SERVICE MENU?** is displayed. Press yes.

Press yes to **PROG. DOWNLOAD?** and yes to **CONFIRM!**

When the message **POWER DOWN/UP** appears, turn the analyzer off, then turn back on after 3 seconds.

The analyzer will display **LOAD NEW PROGRAM - CONFIRM?** Press yes. **SEND <X.KEY> NOW** is displayed. The analyzer is ready to receive data.

- ▲ At this time, you have approximately one minute to begin the download. If a timeout occurs, the analyzer returns to the **LOAD NEW PROGRAM** message. By answering no to this and yes to **GO BACK TO OLD PROGRAM**, the user can discontinue the download process.

Select Send Text file under transfer menu from the tool bar in Hyper Terminal windows.


Browse to the director which stores \*.key and .hex files. Select the .key file, then click at **OPEN**. The analyzer displays Erasing Flash... and then displays **SEND <.HEX> FILE**. Go to Send Text file again. Select Total.hex file, then click at **OPEN**. The analyzer should display a **RECEIVING FILE...** message. This process will take approximately 7-12 minutes.

- ▲ Responses are required during file transfers.


When **DOWNLOAD FINISH** is displayed, press yes. The display should show **CRC CORRECT?** Confirm the number printed is the same as the name of the hex file and answer yes. If the numbers do not agree, answer no. Restart program and repeat the download steps. If numbers still do not match, contact your Bayer representative.


When all steps are complete, the analyzer will reset itself. Upon restarting, the analyzer is running with new revision software.














SERIAL DATA PORT #2










BUBBLE DETECTOR & HEATER #1





**BUTTON READER & PREAMP CONNECTOR**





PWM HEATING & EXT. RAM

# Contacting Bayer Diagnostics

For technical assistance, contact your local authorized representative. For customer service or additional information, contact your local authorized distributor.

## **Addresses and Communication Numbers**

Manufactured for Bayer



Bayer Corporation  
East Walpole, MA 02032-1597



Bayer Diagnostics Mfg (Sudbury) Ltd  
Sudbury, CO10 2XQ, UK

## **Manufactured for Bayer**

Bayer Corporation  
333 Coney Street  
East Walpole, MA 02032-1597  
USA  
508 668-5000

Bayer S.A.  
División Diagnósticos  
Ricardo Gutiérrez 3652  
B1605EHD Munro - Buenos Aires  
Argentina  
+54 (0) 11-4762-7000

Bayer Australia Limited  
ABN 22 000 138 714  
Diagnostics Business Group  
2 Keith Campbell Court  
Scoresby Victoria 3179  
Australia  
+61 (0) 3-9212-8444

Bayer Austria GesmbH  
Geschäftsbereich Diagnostika  
Lerchenfelder Gürtel 9-11  
A-1164 Wien, Austria  
+43 (0) 1-711-46-2424

Bayer s.a.-n.v.  
Division Diagnostics  
Avenue Louise 143 Louizalaan  
1050 Bruxelles-Brussel,  
Belgium  
+32 (0) 2-535-66-81

Bayer S.A.  
Produtos Diagnósticos  
Rua Domingos Jorge 1100  
04779-900 - São Paulo - SP  
Brazil  
+55-11-5694-5166

Bayer Inc.  
Bayer Diagnostics  
Healthcare Division  
77 Belfield Road  
Toronto, Ontario  
Canada M9W 1G6  
416-248-0771

Bayer S.R.O.  
Litvinovska 609/3  
CZ- 190 21 Prague 9 -Prosek  
+420 (0) 2-66101463

Bayer S.A.  
División Diagnóstica  
A.A. 80387  
Av. De las Américas No 57-52  
Santafé de Bogotá, D.C. - Colombia  
+571 (09) 423-4199

Bayer A/S  
Nørgaardsvej 32  
DK-2800 Lyngby, Denmark  
+45-45-23-50-00

Bayer OY  
Suomalaistentie 7  
FIN 02270 Espoo, Finland  
+358-9-887-887

Bayer Diagnostics  
Tour Horizon  
52, quai de Dion Bouton  
92807 Puteaux Cedex, France  
+33 (0) 1-49-06-56-00

Bayer Vital GmbH  
Geschäftsbereich Diagnostics  
Siemensstrasse 3  
D-35463 Fernwald, Germany  
+49 (0) 641- 40-03-170

Bayer Hellas S.A.  
Akakion 54A  
Marousi  
Athens 151 25, Greece  
+30 (0) 1-6883648

Bayer Diagnostics Limited  
20/F Gee Chang Hong Centre  
65 Wong Chuk Hang Road  
Hong Kong  
+852-28147337

Bayer Hungaria Kft.  
H-1012 Budapest, Hungary  
Pálya u.4-6  
+36 (06) 1-212-1540

Bayer Diagnostics India Limited  
589, Sayajipura  
Ajwa Road  
Baroda - 390 019  
Gujarat, India  
+91 (0) 26-5562720

Bayer S.p.A.  
Divisione Diagnostici  
Via Grosio 10/4  
20151 Milano, Italia  
+39-2-3978-1

Bayer Medical Ltd.  
1-19-15, Ebisu  
Shibuya-Ku  
Tokyo 150-0013, Japan  
+81 (0) 3-440-2411

Bayer Diagnostics Korea Limited  
Kye Myung Bldg. 4F  
Myungil-Don, 48-7  
Kangdong-Ku  
Seoul 134070  
Korea  
+82 (0) 2-428-5987

Bayer (Malaysia) Sdn. Bhd  
19th & 20th Floor Wisma MPSA  
P.O. Box 7252  
40708 Shah Alam  
Selangor Darul Ehsan  
Malaysia  
+603 (0) 551-02818

Bayer de México, S.A. de C.V.  
División Diagnósticos  
Via Morelos 330-E  
Santa Clara  
55540 Ecatepec, Estado de México  
México  
Centro de Atención Telefónica (CAT)  
+55 57 28 33 12  
+55 57 28 33 06

Bayer B.V.  
Health Care  
Division Diagnostics  
Energieweg 1  
3641 RT Mijdrecht  
The Netherlands  
+31 (0) 297-280666

Bayer New Zealand Ltd  
Diagnostics Business Group  
3 Argus Place,  
Glenfield,  
Auckland, New Zealand  
+64-800-724-269



Bayer As  
Brennaveien 18  
N-1483 Skytta, Norway  
+47-67-06-86-00

Bayer Sp. Z o.o.  
Al. Jerozolimskie 158  
02-326 Warszawa, Polska  
+48 (0) 22-572-3500

Bayer Diagnostics Europe Ltd  
(sucursal em Portugal)  
Rua Quinta do Pinheiro, 5  
2795-653 Carnaxide, Portugal  
+351-21-416-4227

Bayer Puerto Rico Inc.  
Diagnostics Division  
Victoria Industrial Park  
Building 1  
Carolina, Puerto Rico  
787-752-8989

Bayer (South East Asia) Pte Ltd.  
No. 9 Benoi Sector  
Singapore 629844  
+65-261-3389

Bayer (Pty) Ltd.  
Healthcare Division  
27 Wrench Road  
Isando 1600, South Africa  
+27 (0) 11-921-5048

Química Farmacéutica Bayer, S.A.  
División Diagnósticos  
Calabria, 268  
08029 Barcelona, España  
+34-93-4956500

Bayer AB  
Drakegatan 1  
S-402 24 Göteborg, Sweden  
+46-31-83-98-00

Bayer (Schweiz) AG  
GB DS  
Grubenstrasse 6  
CH-8045 Zürich  
+41 (0) 1-465-81-11

Bayer Taiwan Co. Ltd.  
Diagnostics Division  
11/F 237 Sung Chiang Road  
Taipei 104  
Taiwan  
+886-2-2-5039123

Bayer plc  
Diagnostics Division  
Bayer House  
Strawberry Hill  
Newbury, RG14 1JA  
United Kingdom  
+44 (0) 1635-563000

Bayer Corporation  
511 Benedict Avenue  
Tarrytown, NY 10591-5097  
USA  
914-631-8000

