

## 1.5. Host Query

With the Host Query option enabled, when an unknown patient bar-code is read (patient data not present in the software) a request is sent to the Host computer to send the patient data and related analyses for that unknown code. Moreover the analyzer can send to the Host computer a message containing the sending patients outcome, which may be positive or can contain errors.

The latter can be disabled in the Setup ([see chap. 9, par. 2.2](#))

### Request for patient data + analyses data:

To the host is sent a request as follows:

**HEADER**  
**INFORMATION**  
**TERMINATOR**

Where Header and Terminator are the same used in [1.4. Receive report](#), while Information has the following format.

#### Information [Q]

Information request related to a specific code

Rif. ASTM E1394-91	Field Number	Field Name	Max Len	Content	Note
12.1.1	1	Record Type ID	1	Q	Identifies the message Information (fix)
12.1.5	3	Patient ID	30	Patient code	It is the patient's unique ID code
12.1.13	13	Type	1	O	Patient data request + test

As answer, the LIS will send a pack containing the patient data and the test to be performed for that code, as specified in [1.3. Analyzer Programming](#)

### Sending errors/status to LIS:

A request as follows is sent to the host:

**HEADER**  
**TERMINATOR**

Where Header and Terminator are the same used in [1.4. Receive report](#). Only the Terminator Information field is modified (Rif. 1394-91 13.1.3) as it may contain one of the following codes:

Code	Type
F	Sending result OK
T	The sender has aborted the communication
R	The receiver has aborted the communication
Q	Last pack had an error in contents (usually an error in the analyses codes – too many analyses or unknown analysis code)

**ATTENTION:** when the query mode is activated, make sure that the host is connected, as the operative program has a 60 seconds timeout for each single request. In case the host is not connected, the analyzer might stop working for a long time.