

- If more Order packs are sent for a single patient and only one contains as Test priority the identification as urgent, the whole patient (with all the related analyses) will be placed in the STATs list.

1.4. Receive Report

Header [H]

This message contains only the identification data of the message sender.

Rif. ASTM E1394-91	Field Number	Field Name	Length	Content	Note
7.1.1	1	Record Type ID	1	H	Identifies the message Header (fix)
7.1.2	2	Delimiter definition	4	^&	Delimiters
7.1.5	5	Sender ID	Max 50	Program ID^Version	Analyzer name and software version
7.1.14	14	Date Time	8+2	YYYYMMDD^Time Zone	Sending date + time zone if different from 0

Patient [P]

Sends a single patient basic information

Rif. ASTM E1394-91	Field Number	Field Name	Max Length	Content	Note
8.1.1	1	Record Type ID	1	P	Identifies the message Patient (fix)
8.1.3	3	Patient ID	30	Patient code	It is the patient's unique ID code
8.1.6	6	Patient Name	60	Surname^Name	Patient's surname and name
8.1.8	8	Birth Date	8	YYYYMMDD	Patient's date of birth
8.1.9	9	Patient Sex	3	M, F, C, G1, G2, G3, G4, G5, G6, G7, G8, G9, G10, G11, G12, G13, G14, G15	Patient's sex (group) (Groups from G1 to G15 are the user-defined groups)
8.1.23	22	Note (line 1)	40	Patient's note	Lab notes
8.1.23	23	Note (line 2)	40	Patient's note	Lab notes
8.1.24	24	Admission Date	8	YYYYMMDD	Admission date

Result [R]

Sends a single test result

Rif. ASTM E1394-91	Field Number	Field Name	Max Length	Content	Note
10.1.1	1	Record Type ID	1	R	Identifies the message Result (fix)
10.1.3	3	Test Code	4	Analysis code	Analysis unique identification code
10.1.4	4	Result	10	Result	Result expressed in concentration
10.1.5	5	Units	30	Units	Units
10.1.6	6	Range	20	Min^Max	Normal values range
10.1.7	7	Flag	10	Flags	Analysis flags separated by ^
10.1.9	9	Status	1	F	Result status (fix)
10.1.11	11	Operator	30	Operator	Operator name
10.1.12	12	Start Date	8+2	YYYYMMDD^Time Zone	Test start date + time zone if different from 0

Terminator [L]

This message terminates the communication

Rif. ASTM E1394-91	Field Number	Field Name	Max Length	Content	Note
13.1.1	1	Record Type ID	1	L	Identifies the message Terminator (fix)
13.1.3	3	Terminator Code	1	N	Closes the data pack (fix)