

## 1.6. Report Request

The analyzer does not automatically transmit, but always waits for a specific request from LIS.

### Reports request:

A request as follows is sent to the analyzer:

**HEADER**  
**INFORMATION**  
**TERMINATOR**

Where Header and Terminator are the same used in [1.3. Analyzer Programming](#), while Information has the following format.

### Information [Q]

Report request for a single code or a range of codes

Rif. ASTM E1394-91	Field Number	Field Name	Max Length	Content	Note
12.1.1	1	Record Type ID	1	Q	Identifies the message Information (fix)
12.1.5	5	Universal ID	15+15	Start code^EndCode	Code/codes that identify the patient/patients (see table 7.1)
12.1.13	13	Type	1	F	Request for the final results (fix)

Tabella 7.1

Start Code	End Code	Note
ALL_PATIENTS	-	Request for all available reports
NEXT_PATIENT	-	Request for report of the next patient to be sent to LIS
Patient ID	-	Request for the patient report with patient ID equal to the sent patient ID
Patient ID Start	Patient ID Stop	Request for all patients reports with ID included between Patient ID Start e Patient ID Stop (extremes included)