

## IMPORTANT NOTICE & WARNING

### **Important Notice:**

Any modification to the internal serial protocol is restricted to qualified personnel only. Biotechnica Instruments warrants the correct operation of the internal serial protocol, but does not assume any responsibility for any changes made by any third party. The responsibility for any malfunction arising out of any modifications to the scripts of the internal serial protocol rests with the customer.

### **WARNING**

This information regards the setting up of the barcode for sample tubes identification. The reading of the sample barcode label has the same progression as patient code.

For example: Once a patient code of 15 characters has been entered, then a code of 8 characters followed by 7 empty spaces to reach the 15 characters is sent.

The code read on the barcode label must have the same sequence 8 + 7 for correct detection.