## Communication format

```
Cable: Din9 Standard Serial Cross Cable
      (TX 2 <-> RX 3 Pin Cross, 5 Pin Ground)
      * Another pins must be disconnected.
Protocol: Serial
    -Baud Rate: 115200 baud (fixed)
    -Data Bits: 8-bit
                                 -Stop Bits: 1 -Parity: Disabled
    -Transmission Mode: ASCII
Date Format: 28byte * 18 Line
       1. Start/End Marker : "~" or '0x7e' :1 byte
       2. Each line format (28 byte: include 0x0d, 0x0a)
           (1). 1 line: mm/dd/yy + "" + hh:mm:ss + "" + 0x0d + 0x0a
           (2). 2 line: "Name:"(5byte)+ ""(13 byte) + "Sex:"(4byte)
                      + " "(4 byte) + 0x0d + 0x0a
           (3). 3 line: "Ward:" (5byte) + ""(13 byte) + "Age:" (4byte)
                       + ""(4 byte) + 0x0d + 0x0a
           (4). 4 line: "Strip Name" (10byte) + "" (1byte) + "" (8byte)
                       + 0x0d + 0x0a
           (5). 5-14 line: "URO "(3 byte) +" "(4byte)+ "norm "(4 byte)
                       + " "(3byte)+ " 0.1mg/dl "(12 bytes) + 0x0d + 0x0a
           (6). 15 line: "(SN=0000001)"(12byte) + " "(14byte) + 0x0d+0x0a
           (7). 16 line: "ID(
                                       ) (001)"(23byte)+" "(3byte)
                        + 0x0d + 0x0a
           (8). 17 line: "OP(00000000) LOT(000000) "(25byte) + "~"(1byte)
```

+ 0x0d + 0x0a

```
~DEC/26/2006 11:01:34
```

Name: Sex:

Ward: Age:

Strip Name (10) (0071)

URO norm 0.1mg/dl

GLU - neg

BIL - neg

KET – neg

S.G 1.025

BLD - neg

pH 6.0

PRO - neg

NIT - neg

LEU - neg

(SN=60153025)

ID(

OP(00000000) LOT(000000)  $\sim$